









	Artikel	Specificatie	No.
MLRP	ML RP test set	A ML RP System in an oculus ring	72901
			

MLFILTER Tests		ML Filter Test Set	Choose your own box. Pick a minimum of 4 different ML Filter standard products and the box will be free of charge	
		Empty box		94040
		Flirp test set 1	A flirp stand and the following 3 flirps: - 450PM/511PM, - 527PM/550PM, - Pol1PM/Pol3PM	39920
		Flirp test set 2	A flirp stand and the following 2 flirps: - C1PM/500PM - 400PM/585 CR39	39921
		Flirp test set 3	A flirp stand and the following 5 flirps: - C1PM/500PM - 400PM/585 CR39 - 450PM/511PM - 527PM/550PM - Pol1PM/Pol3PM	39925
		ML Filter Flirp 450PM/511PM	Flirp with ML Filter 450PM & 511PM	39911
		ML Filter Flirp 527PM/550PM	Flirp with ML Filter 527PM & 550PM	39912
		ML Filter Flirp Pol1PM/Pol3PM	Flirp with Pol1PM & Pol3PM	39913
		ML Filter Flirp 400PM/450PM	Flirp with ML Filter 400PM & 450PM	39914
		ML Filter Flirp C1PM/500PM	Flirp with ML Filter C1PM & 500PM	39917
		ML Filter Flirp 400PM/585CR39	Flirp with ML Filter 400PM & 585CR39	39919
		Flirp stand	A table stand for up to 5 flirps	92000
		ML Filter Showframe Box	A compact box with room for up to ten different show frames (based on the sun-two frame) which enables easy testing. The test box contains the following 9 filters in PM material: Pol1, Pol3, C1, 400, 450, 500, 511, 527, 550	39930
			Showframe 585	A showframe with 585 filter in CR39. Not included in the standard Showframebox, but there is a free slot in the box to put it in.
		ML Filter Biocover Test Set	Choose your own box. Pick a minimum of 3 different ML Filter Biocover standard products and the box will be free of charge	
			Empty box	
		Filter in oculusing PM(1)	400, C1, 450, 500, 511, 527 or 550 nm	380**
		Filter in oculusing	585 nm	38058
		Filter in oculusing polarized1	400, C1, 450 (PM), 500, 511 (PM), 527 (PM) or 550 nm	381**
		Filter in oculusing polarized3	400, C1, 450, 500, 511, 527, 550 or 585 nm	383**
Filter in oculusing polarized1		no filter	38100	
Filter in oculusing polarized3		no filter	38300	

** : Initial two figures of the filter colour

(400 = 40 / 450 = 45 / 500 = 50 / 511 = 51 / 527 = 52 / 550 = 55 / 585 = 58)

(1): All PM lenses incl. hard coat